



Immedia BedString

SUPPORTS SITTING UP IN BED

Immedia BedString with handles is mainly intended for the active user to support moving to a sitting position in bed. In the sitting position, the user may use it to lean themselves forward to allow easier application of a sling, positioning of pillows or for dressing. BedString has 8 handles, alternating right and left side of the string to enable the user to find a good grip.



Immedia BedString with 8 handles on alternating sides of the string.

Immedia HandyGrip

ADJUSTABLE GRIP TO AID TURNING IN BED

Immedia HandyGrip provides help with independent turning in bed. It can also be used when self transferring from wheelchair to bed. Soft, compact and easy to hide away. The adjustable buckle allows the length and position of HandyGrip to be varied, providing the user with a good grip in the required position.



Soft, compact and easy to use.

DESCRIPTION	WIDTH X LENGTH		ITEM NO.	SPECIFICATIONS
Immedia	CM	INCH		
~ BedString with handles	8,5 x 170	3 x 66.9	IM429	BEDSTRING: Polyester
~ HandyGrip	5 x 160	2 x 63	IM461	HANDYGRIP: Polyester BELT BUCKLE: Acetal and nylon HANDYGRIP: Polyester
WASHING INSTRUCTIONS				
Do not use fabric softener or bleach.				